


Product datasheet for **GA108641**

Human NIPSNAP3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NIPSNAP3A CRISPRa kit - CRISPR gene activation of human nipsnap homolog 3A
Symbol:	NIPSNAP3A
Components:	GA108641G1, NIPSNAP3A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108641G2, NIPSNAP3A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108641G3, NIPSNAP3A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UFN0</u>
Synonyms:	HSPC299; NIPSNAP4; TASSC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)