

Product datasheet for **GA108648**

Human mSin3A (SIN3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SIN3A CRISPRa kit - CRISPR gene activation of human SIN3 transcription regulator family member A
Symbol:	SIN3A
UniProt ID:	Q96ST3
Synonyms:	WITKOS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)