

Product datasheet for **GA108650**

Human Nectin 3 (NECTIN3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NECTIN3 CRISPRa kit - CRISPR gene activation of human nectin cell adhesion molecule 3
Symbol:	NECTIN3
UniProt ID:	<u>Q9NQ3</u>
Synonyms:	CD113; CDW113; NECTIN-3; PPR3; PRR3; PVRR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)