

Product datasheet for **GA108656**

Human SEC31B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEC31B CRISPRa kit - CRISPR gene activation of human SEC31 homolog B, COPII coat complex component
Symbol:	SEC31B
Components:	GA108656G1, SEC31B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108656G2, SEC31B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108656G3, SEC31B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NQW1</u>
Synonyms:	SEC31B-1; SEC31L2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)