

Product datasheet for **GA108662**

Human TMEM87A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM87A CRISPRa kit - CRISPR gene activation of human transmembrane protein 87A
Symbol:	TMEM87A
UniProt ID:	<u>Q8NBN3</u>
Synonyms:	DKFZp564G2022
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)