

Product datasheet for **GA108669**

Human PARP7 (TIPARP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIPARP CRISPRa kit - CRISPR gene activation of human TCDD inducible poly(ADP-ribose) polymerase
Symbol:	TIPARP
UniProt ID:	Q7Z3E1
Synonyms:	ARTD14; PARP7; pART14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)