


Product datasheet for **GA108677**

Human Tsukushin (TSKU) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSKU CRISPRa kit - CRISPR gene activation of human tsukushi, small leucine rich proteoglycan
Symbol:	TSKU
Components:	GA108677G1, Tsukushin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108677G2, Tsukushin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108677G3, Tsukushin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WUA8
Synonyms:	E2IG4; LRRC54; TSK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)