

Product datasheet for **GA108678**

Human HINFP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HINFP CRISPRa kit - CRISPR gene activation of human histone H4 transcription factor
Symbol:	HINFP
UniProt ID:	<u>Q9BQA5</u>
Synonyms:	HiNF-P; MIZF; ZNF743
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)