

Product datasheet for **GA108680**

Human SNED1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNED1 CRISPRa kit - CRISPR gene activation of human sushi, nidogen and EGF like domains 1
Symbol:	SNED1
Components:	GA108680G1, SNED1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108680G2, SNED1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108680G3, SNED1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TER0</u>
Synonyms:	Snep; SST3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)