

Product datasheet for **GA108684**

Human CLIP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CLIP3 CRISPRa kit - CRISPR gene activation of human CAP-Gly domain containing linker protein 3
Symbol:	CLIP3
UniProt ID:	Q96DZ5
Synonyms:	CLIPR-59; CLIPR59; RSNL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)