

Product datasheet for **GA108686**

Human RNF167 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RNF167 CRISPRa kit - CRISPR gene activation of human ring finger protein 167
Symbol:	RNF167
Components:	GA108686G1, RNF167 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108686G2, RNF167 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108686G3, RNF167 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6Y7</u>
Synonyms:	5730408C10Rik; LP2254; RING105
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)