

Product datasheet for **GA108689**

Human C2CD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C2CD3 CRISPRa kit - CRISPR gene activation of human C2 domain containing 3 centriole elongation regulator
Symbol:	C2CD3
Components:	GA108689G1, C2CD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108689G2, C2CD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108689G3, C2CD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4AC94</u>
Synonyms:	OFD14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)