

Product datasheet for **GA108690**

Human DAK (TKFC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TKFC CRISPRa kit - CRISPR gene activation of human triokinase and FMN cyclase
Symbol:	TKFC
UniProt ID:	Q3LXA3
Synonyms:	DAK; NET45
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)