

## Product datasheet for **GA108692**

### Human SPATS2L activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | SPATS2L CRISPRa kit - CRISPR gene activation of human spermatogenesis associated serine rich 2 like                                                                                                                                                                                                              |
| Symbol:       | SPATS2L                                                                                                                                                                                                                                                                                                          |
| Components:   | GA108692G1, SPATS2L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108692G2, SPATS2L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108692G3, SPATS2L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9NUQ6</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | DNATP6; SGNP                                                                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)