

Product datasheet for **GA108693**

Human Teneurin 4 (TENM4) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | TENM4 CRISPRa kit - CRISPR gene activation of human teneurin transmembrane protein 4 |
| Symbol: | TENM4 |
| UniProt ID: | Q6N022 |
| Synonyms: | Doc4; ODZ4; Ten-M4; TNM4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)