

Product datasheet for **GA108698**

Human LRIG1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRIG1 CRISPRa kit - CRISPR gene activation of human leucine rich repeats and immunoglobulin like domains 1
Symbol:	LRIG1
UniProt ID:	Q96JA1
Synonyms:	LIG-1; LIG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)