

## Product datasheet for **GA108705**

### Human PLEKHG3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHG3 CRISPRa kit - CRISPR gene activation of human pleckstrin homology and RhoGEF domain containing G3
Symbol:	PLEKHG3
UniProt ID:	<a href="#">A1L390</a>
Synonyms:	ARHGEF43; KIAA0599
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

