


## Product datasheet for **GA108714**

### Human SYT homolog 1 (SS18L1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SS18L1 CRISPRa kit - CRISPR gene activation of human SS18L1 subunit of BAF chromatin remodeling complex
Symbol:	SS18L1
Components:	GA108714G1, SYT homolog 1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108714G2, SYT homolog 1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108714G3, SYT homolog 1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q75177</a></u>
Synonyms:	CREST; LP2261
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)