

Product datasheet for **GA108724**

Human Dynamin 3 (DNM3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNM3 CRISPRa kit - CRISPR gene activation of human dynamin 3
Symbol:	DNM3
UniProt ID:	Q9UQ16
Synonyms:	Dyna III
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)