

Product datasheet for **GA108735**

Human LSM14A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LSM14A CRISPRa kit - CRISPR gene activation of human LSM14A mRNA processing body assembly factor
Symbol:	LSM14A
Components:	GA108735G1, LSM14A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108735G2, LSM14A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108735G3, LSM14A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8ND56
Synonyms:	C19orf13; FAM61A; RAP55; RAP55A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)