

Product datasheet for **GA108736**

Human RTL8A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RTL8A CRISPRa kit - CRISPR gene activation of human retrotransposon Gag like 8A
Symbol:	RTL8A
Components:	GA108736G1, RTL8A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108736G2, RTL8A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108736G3, RTL8A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BWD3</u>
Synonyms:	CXX1b; MAR8A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)