

Product datasheet for **GA108737**

Human POLDIP2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	POLDIP2 CRISPRa kit - CRISPR gene activation of human DNA polymerase delta interacting protein 2
Symbol:	POLDIP2
UniProt ID:	Q9Y2S7
Synonyms:	p38; PDIP38; POLD4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)