


## Product datasheet for **GA108748**

### Human WDR21A (DCAF4) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | DCAF4 CRISPRa kit - CRISPR gene activation of human DDB1 and CUL4 associated factor 4                                                                                                                                                                                                                         |
| Symbol:       | DCAF4                                                                                                                                                                                                                                                                                                         |
| Components:   | GA108748G1, WDR21A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108748G2, WDR21A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108748G3, WDR21A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8WV16</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | WDR21; WDR21A                                                                                                                                                                                                                                                                                                 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)