


Product datasheet for **GA108749**

Human PTPN20B (PTPN20) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPN20 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 20
Symbol:	PTPN20
UniProt ID:	Q4JDL3
Synonyms:	bA42B19.1; bA142I17.1; CT126; PTPN20A; PTPN20B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)