

Product datasheet for **GA108750**

Human C1orf77 (CHTOP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CHTOP CRISPRa kit - CRISPR gene activation of human chromatin target of PRMT1
Symbol:	CHTOP
UniProt ID:	Q9Y3Y2
Synonyms:	C1orf77; FL-SRAG; FOP; pp7704; SRAG; SRAG-3; SRAG-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)