

Product datasheet for **GA108759**

Human ARH (LDLRAP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LDLRAP1 CRISPRa kit - CRISPR gene activation of human low density lipoprotein receptor adaptor protein 1
Symbol:	LDLRAP1
UniProt ID:	Q5SW96
Synonyms:	ARH; ARH1; ARH2; FHCB1; FHCB2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)