

## Product datasheet for **GA108766**

### Human TRPC4AP activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TRPC4AP CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily C member 4 associated protein
Symbol:	TRPC4AP
UniProt ID:	<a href="#">Q8TEL6</a>
Synonyms:	C20orf188; PPP1R158; TRRP4AP; TRUSS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)