


Product datasheet for **GA108768**

Human Testin (TES) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TES CRISPRa kit - CRISPR gene activation of human testin LIM domain protein
Symbol:	TES
Components:	GA108768G1, Testin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108768G2, Testin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108768G3, Testin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UGI8</u>
Synonyms:	TESS; TESS-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)