


Product datasheet for **GA108777**

Human ZNF337 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF337 CRISPRa kit - CRISPR gene activation of human zinc finger protein 337
Symbol:	ZNF337
Components:	GA108777G1, ZNF337 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108777G2, ZNF337 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108777G3, ZNF337 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y3M9</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)