


Product datasheet for **GA108790**

Human TMEM251 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM251 CRISPRa kit - CRISPR gene activation of human transmembrane protein 251
Symbol:	TMEM251
Components:	GA108790G1, TMEM251 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108790G2, TMEM251 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108790G3, TMEM251 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N6I4</u>
Synonyms:	C14orf109
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)