

Product datasheet for **GA108794**

Human PTPN22 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | PTPN22 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 22 |
| Symbol: | PTPN22 |
| UniProt ID: | Q9Y2R2 |
| Synonyms: | LYP; LYP1; LYP2; PEP; PTPN8 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)