

Product datasheet for **GA108803**

Human PHGDH activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHGDH CRISPRa kit - CRISPR gene activation of human phosphoglycerate dehydrogenase
Symbol:	PHGDH
UniProt ID:	<u>O43175</u>
Synonyms:	3-PGDH; 3PGDH; HEL-S-113; NLS; NLS1; PDG; PGAD; PGD; PGDH; PHGDHD; SERA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)