

Product datasheet for **GA108817**

Human KLHL3 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | KLHL3 CRISPRa kit - CRISPR gene activation of human kelch like family member 3 |
| Symbol: | KLHL3 |
| Components: | GA108817G1, KLHL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108817G2, KLHL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108817G3, KLHL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q9UH77</u> |
| Synonyms: | PHA2D |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)