

Product datasheet for **GA108818**

Human **KCN G2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCN G2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily G member 2
Symbol:	KCN G2
UniProt ID:	Q9UJ96
Synonyms:	KCNF2; KV6.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)