


Product datasheet for **GA108838**

Human VPS33B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	VPS33B CRISPRa kit - CRISPR gene activation of human VPS33B late endosome and lysosome associated
Symbol:	VPS33B
Components:	GA108838G1, VPS33B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108838G2, VPS33B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108838G3, VPS33B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H267</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)