

Product datasheet for **GA108839**

Human Tin2 (TINF2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TINF2 CRISPRa kit - CRISPR gene activation of human TERF1 interacting nuclear factor 2
Symbol:	TINF2
UniProt ID:	Q9BSI4
Synonyms:	DKCA3; TIN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)