

## Product datasheet for **GA108864**

### Human PTPN18 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | PTPN18 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 18 |
| Symbol:       | PTPN18   |
| UniProt ID:   | <a href="#">Q99952</a>   |
| Synonyms:     | BDP1; PTP-HSCF   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                |

#### Product images:



[View online >](#)