


Product datasheet for **GA108891**

Human **IL36RN** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	IL36RN CRISPRa kit - CRISPR gene activation of human interleukin 36 receptor antagonist
Symbol:	IL36RN
UniProt ID:	<u>Q9UBH0</u>
Synonyms:	FIL1; FIL1(DELTA); FIL1D; IL-36Ra; IL1F5; IL1HY1; IL1L1; IL1RP3; IL36RA; PSORP; PSORS14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)