

Product datasheet for **GA108905**

Human **STK23 (SRPK3)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SRPK3 CRISPRa kit - CRISPR gene activation of human SRSF protein kinase 3
Symbol:	SRPK3
UniProt ID:	<u>Q9UPE1</u>
Synonyms:	MSSK-1; MSSK1; STK23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)