

Product datasheet for **GA108919**

Human ELP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ELP4 CRISPRa kit - CRISPR gene activation of human elongator acetyltransferase complex subunit 4
Symbol:	ELP4
Components:	GA108919G1, ELP4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108919G2, ELP4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108919G3, ELP4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96EB1</u>
Synonyms:	AN; C11orf19; dj68P15A.1; hELP4; PAX6NEB; PAXNEB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)