

Product datasheet for **GA108941**

Human SH3YL1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3YL1 CRISPRa kit - CRISPR gene activation of human SH3 and SYLF domain containing 1
Symbol:	SH3YL1
Components:	GA108941G1, SH3YL1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108941G2, SH3YL1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108941G3, SH3YL1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96HL8</u>
Synonyms:	RAY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)