

Product datasheet for **GA108945**

Human SMR3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMR3A CRISPRa kit - CRISPR gene activation of human submaxillary gland androgen regulated protein 3A
Symbol:	SMR3A
UniProt ID:	Q99954
Synonyms:	P-B1; PBI; PRL5; PROL5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)