

Product datasheet for **GA108965**

Human **KCNV1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNV1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily V member 1
Symbol:	KCNV1
UniProt ID:	Q6PIU1
Synonyms:	HNKA; KCNB3; KV2.3; KV8.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)