

## Product datasheet for **GA108966**

### Human C2orf24 (CNPPD1) activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| <b>Description:</b>  | CNPPD1 CRISPRa kit - CRISPR gene activation of human cyclin Pas1/PHO80 domain containing 1                                                                                                                                                                                                                       |
| <b>Symbol:</b>       | CNPPD1                                                                                                                                                                                                                                                                                                           |
| <b>Components:</b>   | GA108966G1, C2orf24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108966G2, C2orf24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108966G3, C2orf24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q9BV87</a></u>                                                                                                                                                                                                                                                                                    |
| <b>Synonyms:</b>     | C2orf24; CGI-57                                                                                                                                                                                                                                                                                                  |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:

