

Product datasheet for **GA108980**

Human PKD2L2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKD2L2 CRISPRa kit - CRISPR gene activation of human polycystin 2 like 2, transient receptor potential cation channel
Symbol:	PKD2L2
UniProt ID:	Q9NZM6
Synonyms:	TRPP5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)