

Product datasheet for **GA108992**

Human LAMP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LAMP3 CRISPRa kit - CRISPR gene activation of human lysosomal associated membrane protein 3
Symbol:	LAMP3
UniProt ID:	Q9UQV4
Synonyms:	CD208; DC-LAMP; DC LAMP; DCLAMP; LAMP; LAMP-3; TSC403
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)