

Product datasheet for **GA108999**

Human CD57 (B3GAT1) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | B3GAT1 CRISPRa kit - CRISPR gene activation of human beta-1,3-glucuronyltransferase 1 |
| Symbol: | B3GAT1 |
| UniProt ID: | Q9P2W7 |
| Synonyms: | CD57; GLCATP; GLCUATP; HNK1; LEU7; NK-1; NK1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)