

Product datasheet for **GA109004**

Human **KCNMB3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNMB3 CRISPRa kit - CRISPR gene activation of human potassium calcium-activated channel subfamily M regulatory beta subunit 3
Symbol:	KCNMB3
UniProt ID:	Q9NPA1
Synonyms:	BKBETA3; HBETA3; K(VCA)BETA-3; KCNMB2; KCNMBL; SLO-BETA-3; SLOBETA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)