

Product datasheet for **GA109012**

Human DMAC2L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DMAC2L CRISPRa kit - CRISPR gene activation of human distal membrane arm assembly complex 2 like
Symbol:	DMAC2L
Components:	GA109012G1, DMAC2L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109012G2, DMAC2L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109012G3, DMAC2L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99766</u>
Synonyms:	ATPW; FB; HSU79253
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)