

Product datasheet for **GA109029**

Human **KCNH5** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNH5 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 5
Symbol:	KCNH5
UniProt ID:	<u>Q8NCM2</u>
Synonyms:	EAG2; H-EAG2; hEAG2; Kv10.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)